


phlebology measuring charts

MEASURING INSTRUCTIONS

1. Select product/style from measuring chart. 2. Measure at indicated locations on diagram. 3. Locate proper product size based on measurements.

mediven® assure 15–20, 20–30, 30–40 mmHg

| <i>Circumference</i> | <i>SMALL</i> | <i>MEDIUM</i> | <i>LARGE</i> | <i>XLARGE</i> | <i>XXLARGE</i> |
|---------------------------------|---|-----------------------|----------------------|-----------------------|-----------------------|
| 3 <i>PANTY</i> | 20"–26" 51–66 cm | 22"–28" 56–71 cm | 24"–30" 61–76 cm | 26"–32" 66–81 cm | N/A |
| 2 <i>EXTRA-WIDE CALF</i> | N/A | N/A | 17"–20" 43–51 cm | 18"–22" 46–56 cm | 20"–24" 51–61 cm |
| 2 <i>CALF</i> | 11"–15" 28–38 cm | 12½"–16½" 31–42 cm | 13½"–18" 34–46 cm | 14½"–19½" 37–50 cm | 18"–22" 46–56 cm |
| 1 <i>ANKLE</i> | 7"–8½" 18–21.5 cm | 8½"–9¾" 21.5–25 cm | 9¾"–11" 25–28 cm | 11"–12½" 28–32 cm | 12½"–14½" 32–37 cm |
| <i>LEG LENGTH</i> | C <i>PANTY & THIGH</i> ▶ <i>STANDARD</i> : 28"–34" (71–86.5 cm) <i>PETITE</i> : 23"–28" (58.5–71 cm) | | | | |



For All Panty & Thigh Styles:

Measure **1** **2** **3** and **C**