

SIGVARIS USA CUSTOM ORDERING INSTRUCTIONS



***Please do not type in all capital letters

***All boxes in this section need to be completed so the proper stockings can be ordered for your patient

***Please be sure to include your choice for Grip-top. Available styles are:

- Sensinnov Grip-top (5cm width)
- Knobbed Grip-top (5cm width)
- Narrow Knobbed Grip-top (Calf only, 2.5cm width)

***Bodyform Pantyhose have compression from Ankle to Waist

***Pantyhose for men is equipped with Fly in front (choose Left, Right, or Horizontal Fly in Remarks)

***If you need swatches of material and color choices or Body Pens for measuring, please contact Melissa Gill at email or phone below

***Any questions can be directed to
ACOLS
Sales@acols.com
800-863-5935 Option #1

***Choose your compression level below and click the corresponding box under COMFORT, COTTON, or TRADITIONAL Series

SIGVARIS Custom stockings use German Standard RAL-GZG for compression ranges.

1 = 18–21mmHg
 2 = 23–32mmHg
 3 = 34–46mmHg
 4 = at least 49mmHg

		Comfort		Cotton			Traditional		
		1	2	1	2	3	2	3	4
		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Models									
A-D Calf	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-F Half-thigh	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-G Thigh	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-G Thigh with waist attachment	left	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	right	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose Maternity		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose Bodyform		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose for man		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
One leg Pantyhose		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
One leg Pantyhose slip form		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Colors									
Skin		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Caramel		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Grey		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Black		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Savanna		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Nature		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Beige		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Alternatives: for TRADITIONAL compression up to the cG level (without soft-part)
 * This option is ideal for lymphedema patients.

***Choose either Pairs or Pieces for only 1 leg

Pairs	Pieces
<input type="checkbox"/>	<input type="checkbox"/>

***Click on Model of stocking below COTTON, COMFORT, or TRADITIONAL Series

***Click on the Color choice

***TRADITIONAL Series only available for Alternatives

SIGVARIS USA CUSTOM ORDERING INSTRUCTIONS



***All measurements should be done in CENTIMETERS

***Measuring should be done in the morning and patient should be free of edema to get correct measurements. Patient should be accessible from all sides.

***Measurements should be done on bare skin - do not measure over compression stockings or clothes.

***Measuring tape should not be too tight or loose.


***Measurements should be taken from distal to proximal, i.e., starting at foot.

***Recommend taking measurements in a vertical position.

***Length measurements should be taken straight up the leg and not to the contour of the leg.

***Don't forget the Foot Length for Open or Closed Toe.

***Any questions please call or email Kristin Fields
Direct: (800) 863-5935 Option #1
Sales@acols.com



**Measuring Requirements for Custom Stockings

Calf = cA through cD & lB through lD

Thigh = cA through cG & lB through lG

Thigh w/ Waist Attachment = cA through cG & lB through lT - also include cH & cT & IK-T front/back

Any Pantyhose = cA through cG & lB through lT - also include cH & cT & IK-T front/back

SIGVARIS USA CUSTOM ORDER FORM

Dealer:	Account / Customer #: 11700	Remarks (indicate special order alterations):
Patient Name:	Name: ACOLS	
Patient Gender: F ___ M ___	Address 1: 11632 High St Ste A	
Purchase Order #:	Address 2:	
Reorder #:	Apt./Suite/Other:	
Contact:	City: Sebastian	
Contact Phone:	State: FL Zip Code: 32958	
Contact Email:		

SIGVARIS Custom stocking use German Standard RAL-GZG for compression ranges.

- 1 = 18-21mmHg
- 2 = 23-32mmHg
- 3 = 34-46mmHg
- 4 = at least 49mmHg

Comfort	Cotton	Traditional	Pair	Pieces
1 2	1 2 3	2 3 4	<input type="checkbox"/>	<input type="checkbox"/>
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		

Models

A-D Calf	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-F Half-thigh	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-G Thigh	with soft top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	with grip top	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-G Thigh with waist Attachment	left	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	right	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose Maternity		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose Body form		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
A-T Pantyhose for man		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
One leg Pantyhose		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
One leg Pantyhose slip form		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Colours

Skin	<input type="radio"/>
Caramel	<input type="radio"/>
Grey	<input type="radio"/>
Black	<input type="radio"/>
Savanna	<input type="radio"/>
Nature	<input type="radio"/>
Beige	<input type="radio"/>

Alternatives: for TRADITIONAL compression up to the cG level (without soft-part)
 * This option is ideal for lymphedema patients.

Leg measurements

Centimeters only

Measurements must be taken in the morning on a leg free from edema. We cannot be held responsible for incorrect measurements given. Product liability cannot be accepted if the label has been removed.

Circumference Right Leg	Waist	Circumference Left Leg	Lengths
ℓ K-T front		cT	ℓ T
ℓ K-T back	Hip	cH	ℓ H
cG		cG	ℓ G
cF		cF	ℓ F
cE		cE	ℓ E
cD		cD	ℓ D
cC		cC	ℓ C
cB'		cB'	ℓ B'
cB		cB	ℓ B
cY		cY	
cA		cA	head of metatarsus
			<input type="radio"/> open toe
			<input type="radio"/> closed toe

Shoe Size

Grip-tops
<input type="radio"/> Sensinno [®] grip top (approx. 5cm)
<input type="radio"/> Knobbed grip top (approx. 5cm)
<input type="radio"/> Narrow Knobbed grip top (approx. 5cm)