


phlebology measuring charts

MEASURING INSTRUCTIONS

1. Select product/style from measuring chart. 2. Measure at indicated locations on diagram. 3. Locate proper product size based on measurements.

mediven® assure 15–20, 20–30, 30–40 mmHg

<i>Circumference</i>	<i>SMALL</i>	<i>MEDIUM</i>	<i>LARGE</i>	<i>XLARGE</i>	<i>XXLARGE</i>
3 <i>PANTY</i>	20"–26" 51–66 cm	22"–28" 56–71 cm	24"–30" 61–76 cm	26"–32" 66–81 cm	N/A
2 <i>EXTRA-WIDE CALF</i>	N/A	N/A	17"–20" 43–51 cm	18"–22" 46–56 cm	20"–24" 51–61 cm
2 <i>CALF</i>	11"–15" 28–38 cm	12½"–16½" 31–42 cm	13½"–18" 34–46 cm	14½"–19½" 37–50 cm	18"–22" 46–56 cm
1 <i>ANKLE</i>	7"–8½" 18–21.5 cm	8½"–9¾" 21.5–25 cm	9¾"–11" 25–28 cm	11"–12½" 28–32 cm	12½"–14½" 32–37 cm
<i>LEG LENGTH</i>	C <i>PANTY & THIGH</i> ▶ <i>STANDARD</i> : 28"–34" (71–86.5 cm) <i>PETITE</i> : 23"–28" (58.5–71 cm)				



For All Panty & Thigh Styles:

Measure **1** **2** **3** and **C**