


# phlebology measuring charts

## MEASURING INSTRUCTIONS

1. Select product/style from measuring chart. 2. Measure at indicated locations on diagram. 3. Locate proper product size based on measurements.

### mediven® active 15–20, 20–30 mmHg

Circumference	I	II	III	IV	V	VI	VII
<b>2 CALF</b>	11"–14" 28–35.5 cm	11½"–15¼" 30–39 cm	12¾"–16½" 32.5–42 cm	13¾"–17¾" 35–45 cm	14½"–19" 37–48.5 cm	15½"–20¼" 39.5–51.5 cm	16¾"–21½" 41.5–54.5 cm
<b>1 ANKLE</b>	7"–7½" 18–20 cm	8"–8½" 20.5–22.5 cm	9"–9½" 23–25 cm	10"–10½" 25.5–27.5 cm	11"–11½" 28–30 cm	12"–13" 30.5–33 cm	13¾"–14¼" 33.5–36 cm
<b>LEG LENGTH</b>	<b>B CALF</b> ▶ <b>STANDARD:</b> 15"–19" (38–47.5 cm) <b>TALL:</b> >18" (>47.5 cm)						



For All Calf Styles:  
Measure **1** **2** and **B**