



## PRIMA



### Band Size

centimetres

66-71		
76-81	<b>S</b>	<b>M</b>
86-91	<b>M</b>	
96-102	<b>L</b>	
107-112	<b>XL</b>	<b>L</b>
117-122	<b>2XL</b>	
127-132	<b>2XL With Extender</b>	

HUGGER PRIMA Prima EXTENDER

High to Medium Compression Fits your HuggerPRIMA

# FIT CHART

**You just need one measurement :)**

**Really!**

Band size is your current bra size (North America) or your rib band size shown here ie: 90 cm = M



#### In between sizes? No worries :)

HuggerPRIMA - Upsize to the nearest size.  
ie: 94 cm = size L

HuggerVIDA - Downsize to the nearest size.  
ie: 94 cm = size M

#### DD cup or higher? Sensitive to compression?

For the HuggerPRIMA, we recommend that you upsize one size for comfort and wearability or add an EXTENDER. With the 360 seamless knit mapped compression, you can still maintain good compression by upsizing and, most importantly, gain patient compliance with comfortability.

#### Use the EXTENDER for flexibility :)

Our patent-pending Extenders are made to easily zip into the Hugger to vary the compression & breathability instantly • accommodate fluctuating cup/body size or additional compression pads/bandages • allow for 24/7 post-surgical & lymphedema comfort • allow you to easily adjust the size and fit of the Hugger.



## VIDA



### Band Size

centimetres

66-71	<b>XS</b>	
76-81	<b>S</b>	<b>M</b>
86-91	<b>M</b>	
96-102	<b>L</b>	
107-112	<b>XL</b>	<b>L</b>
117-122	<b>2XL</b>	
127-132	<b>2XL With Extender</b>	

HUGGER VIDA VIDA EXTENDER

Medium to Low Compression Fits your HuggerVIDA

# FIT CHART

**You just need one measurement :)**

**Really!**

Band size is your current bra size (North America) or your rib band size shown here ie: 90 cm = M



#### In between sizes? No worries :)

HuggerPRIMA - Upsize to the nearest size.  
ie: 94 cm = size L

HuggerVIDA - Downsize to the nearest size.  
ie: 94 cm = size M

#### DD cup or higher? Sensitive to compression?

For the HuggerPRIMA, we recommend that you upsize one size for comfort and wearability or add an EXTENDER. With the 360 seamless knit mapped compression, you can still maintain good compression by upsizing and, most importantly, gain patient compliance with comfortability.

#### Use the EXTENDER for flexibility :)

Our patent-pending Extenders are made to easily zip into the Hugger to vary the compression & breathability instantly • accommodate fluctuating cup/body size or additional compression pads/bandages • allow for 24/7 post-surgical & lymphedema comfort • allow you to easily adjust the size and fit of the Hugger.