


# circaid reduction kit sizing information for leg



| circaid reduction kit - upper leg |         |          |
|-----------------------------------|---------|----------|
| Circumference in cm               | regular | wide     |
| cG                                | < 80    | < 125    |
| cE2                               | < 65    | < 110    |
| Length in cm                      | short   | standard |
| l D1 - G                          | < 37    | 37 - 41  |
|                                   |         | long     |
|                                   |         | > 41     |

## Each upper leg kit includes:

- 1 reduction upper leg component
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

| circaid reduction kit - lower leg |         |          |
|-----------------------------------|---------|----------|
| Circumference in cm               | regular | wide     |
| cC                                | < 65    | < 80     |
| cB                                | < 55    | < 65     |
| Length in cm                      | short   | standard |
| l mall - D1                       | < 32    | 32 - 36  |
|                                   |         | long     |
|                                   |         | > 36     |

## Each lower leg kit includes:

- 1 reduction lower leg component
- 1 pair reduction kit undersleeve leg
- 1 reduction shelf strap
- 6 hook and loop stays
- 1 Built-In-Tension system guide card
- 1 paper measuring tape
- 1 direction for use

| circaid reduction kit - knee |               |
|------------------------------|---------------|
| Circumference in cm          | regular       |
| cE2                          | < 110         |
| cE                           | < 102         |
| cC                           | < 80          |
| one length available:        | standard 30cm |

## Each knee kit includes:

- 1 reduction knee component
- 2 hook and loop knee spines
- 1 paper measuring tape
- 1 direction for use
- 1 reduction kit shelf strap